

# Online FT-NIR analysis of propylene oxide during polyether polyol production process



Detection of ppm levels of propylene oxide for safe polyether polyol production process.

—  
01 FTPA200-260CH FT-NIR  
multichannel chemical  
process analyzer

## Overview

Polyether polyol is a key component used in the producing of polyurethane; it is synthesized through the reaction of a glycol and an organic oxide. The reaction process often involves a special double metal cyanide catalyst (DMC), which renders the process virtually wasteless and eliminates the use of highly reactive and toxic isocyanate compounds. The most common type of organic oxide used in the DMC catalyst production process of polyether polyol is propylene oxide.

Despite being less toxic and dangerous than isocyanate compounds, propylene oxide still remains an extremely flammable substance with acute toxicity and carcinogenic effects. Thereby, detection of propylene oxide at ppm levels becomes crucial to ensure a safe polyether polyol product at the reactor outlet.

The typical detection method of propylene oxide in polyether polyol product uses gas chromatography, which requires off-line sampling and extraction of the product prior to measurement. Aside from incurring a long detection response time due to the inherent laboratory method, this approach also presents a health risk to personnel.

For such applications, using Fourier transform near-infrared (FT-NIR) spectroscopy with an insertion measuring probe or flow-through cell becomes advantageous because it allows fast and reliable online or inline monitoring of low concentration levels of propylene oxide, thus reducing operators' health risk during the polyol production process.

## Method

Instrument: FTPA2000-260CH multichannel chemical process analyzer

Detector: Room temperature InGaAs 2.6  $\mu\text{m}$

Sampling technique: PSD series 500 single pass transmission probe

Analysis temperature: 150°C


Resolution: 16  $\text{cm}^{-1}$

Number of scans: 64 to 128




01 PSD series 500 single pass transmission probe schematic

02 Propylene oxide concentration calibration results (FT-NIR predicted values against actual values)



## Results



## Calibration results

Property	Range (ppm)	$r^2$	SECV (ppm)
Propylene oxide	0.012 – 1.208	0.986	0.032

## Conclusion

Concentrations of propylene oxide measured with FT-NIR process analyzer are in agreement with the laboratory reference values. This study demonstrates that FT-NIR process analysis can be successfully and reliably used to quantify and monitor in real-time the concentration of propylene oxide in the DMC-catalyzed production process of polyether polyol.

ABB Inc.  
Measurement & Analytics  
3400, rue Pierre-Ardouin, Québec, Québec  
G1P 0B2 Canada

1 800 858-3847 (North America)  
Tel.: +1 418-877-2944 (other countries)  
Fax: +1 418-877-2834  
Email: ftir@ca.abb.com

[abb.com/analytical](http://abb.com/analytical)

We reserve the right to make technical changes or modify the contents of this document without prior notice. With regard to purchase orders, the agreed particulars shall prevail. ABB does not accept any responsibility whatsoever for potential errors or possible lack of information in this document.

We reserve all rights in this document and in the subject matter and illustrations contained therein. Any reproduction, disclosure to third parties or utilization of its contents – in whole or in parts – is forbidden without prior written consent of ABB.  
Copyright© 2024 ABB. All rights reserved