

Malvern Particle Size Analyzer – Insitec Driver/Interface

This document describes the driver/interface for the Malvern Insitec Voyager Pharma process particle size analyzer.


Product Description

The interface provides data collection, instrument parameter control and sampling system control from xPAT. This allows control of sample conditioning and acquisition of background calibration data as well as sample data; e.g. particle size histograms and associated statistics.

The Insitec Voyager Pharma is an integrated analyzer with a built in sample conditioning system and Insitec particle size system. The Insitec Voyager is compatible with Insitec X, D or T optical heads.

The interface is made up of the Malvern RTSizer 7.2.0 software that allows acquisition of particle size distribution histograms, an OPC/UA server that Interfaces to RTSizer, a CANopen I/O interface to control the sampling system, the xPAT analyzer service, the xPAT configuration template for Malvern Insitec and the xPAT object type for Malvern Insitec.

Manual control and status display of the analyzer is provided by a standard faceplate associated with each instance of the analyzer on the xPAT workplace.



Malvern Insitec Voyager Pharma process particle size analyzer




Photo courtesy of Malvern Instruments Ltd.

Specifications

| | |
|---|---|
| Analyzer Class | Particle size analyzer |
| Subclass | Laser Diffraction |
| Interface | OPC/UA for particle size analyzer / CANopen over CANbus for sampling system control |
| Compatibility | Insitac T,D,X,L,S,SX Optical measurement heads Insitac Voyager Pharma integrated system (Sampling system & optical measurement head) |
| Throughput | Max 1 sample every 5 seconds |
| Histogram Size | Nominal 128 points |
| Control Parameters | Read/write access to all parameters |
| Channels | 1 Channel |
| Averaging Time | Averaging time for background / sample measurements in seconds |
| Backflush | Cleaning sequence for the particle size sensor |
| Introduce Sample | Control sampling system to allow a sample into the particle size sensor |
| Signal Processing Parameters | Settings for signal processing of raw data |
| Particle Size Range | Minimum and maximum particle size in microns |
| Media Parameters | Parameters that describe the carrier media; e.g. expected particle density, index of refraction |
| Analysis Parameters | Parameters to interpret the raw scattering data; e.g. multiple scattering on/off |
| Correlation Parameters | Curve fit parameters from scattering data to particle size distribution |
| Faceplate Status Indicators – analyzer | |
| Connection Status | Status of link to analyzer: good or bad |
| Analyzer Status | Status of analyzer: good or bad |
| Faceplate Status Indicators per channel | |
| Acquisition Status | Idle, Sample starting, Sample |
| Data Collection Status | Normal, Maintenance, Fault |
| Faceplate Commands per channel | Collect Sample, Collect Background |
| Control Type | xPAT provides start/stop signal |
| Data Acquisition | Collect dark and light background scattering for calibration and histogram of particle size distribution |
| Data Analysis | 16 statistics computed from particle size distribution |
| Calibration | Collect dark and light background |
| Validation | Operational Qualification (OQ) of analyzer supported, requires measurement of sample with known particle size distribution. Performance Qualification (PQ) implemented by method specific configuration |
| Health Monitoring | Monitor analyzer hardware status; e.g. connections to analyzer and sampling system, monitor analyzer internal parameters |
| Asset Management | Not implemented |

3EUS094897

For more information on the Malvern Particle Size Analyzers please visit www.malvern.com/insitac. For more information on ABB Life Sciences solutions visit www.abb.com/lifesciences.

Note:

We reserve the right to make technical changes or modify the contents of this document without prior notice. With regard to purchase orders, the agreed particulars shall prevail. ABB does not accept any responsibility whatsoever for potential errors or possible lack of information in this document. We reserve all rights in this document and in the subject matter and illustrations contained therein. Any reproduction, disclosure to third parties or utilization of its contents - in whole or in part - is forbidden without prior written consent of ABB. The IndustrialIT wordmark, Aspect Objects, and all above mentioned names in the form xxxxxxIT are registered or pending trademarks of ABB. All rights to other trademarks reside with their respective owners. © Copyright 2009 ABB. All rights reserved.

Power and productivity
for a better world™

