pH BUFFER POWDER INSTRUCTIONS

For many years, Buffer Tablets have been used to standardise pH Meters. The new Buffer Powders represent an important advance in the preparation of standard solutions since they are more accurate and dissolve more easily than tablets. The pure dry chemicals used are sealed into a laminate of polythene and aluminium foil, which keeps out moisture and gases that can cause deterioration of the buffer powder.

1

Before opening, the powder must be shaken down to the bottom of the sachet. This may be done conveniently by holding the sachet upright and tapping the bottom edge on a hard surface, i.e., a bench top. Alternatively, the sachet may be held upright and the top flicked with a finger nail.

2

The corner of the sachet should then be cut off and the contents emptied into the vessel used for dissolving the powder. The sachet should be tapped to loosen any powder adhering to the inside.

3

The water should be freshly distilled or demineralised. Distillate from hard water may be heavily contaminated with carbon dioxide, which should be removed by aeration or boiling prior to use for 9pH buffer. An accuracy of 2% in the volume of water used is sufficient.

4

Neutral salts affect the pH value; undue contamination with salt bridge solution should be avoided.

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Each sachet makes 200 ml of solution.
The pH value is in accordance with NIST, DIN and BS1647 pt. 2, 1984 specifications i.e. 4.01 ±0.02pH at 25°C.
pH values at other temperatures are:

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Each sachet makes 200 ml of solution.
The pH value is 7.00 ±0.04pH at 25°C.
pH values at other temperatures are:

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Each sachet makes 200 ml of solution.
The pH value is in accordance with BS1647 pt. 2, 1984 specification i.e. 9.18 ±0.02pH at 25°C.
pH values at other temperatures are:

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7pH

Di-Sodium Hydrogen Phosphate and Potassium
Dihydrogen Phosphate
(Na₂HPO₄ and NaH₂PO₄)

°C	0	10	20	25	30	40
pH	7.11	7.06	7.01	7.00	6.98	6.97
°C	50	60	70	80	90	
pH	6.97	6.97	6.99	7.03	7.08	

9pH

Sodium Tetraborate
(Na₂B₄O₇)

°C	0	10	20	25	30	40
pH	9.48	9.35	9.23	9.18	9.13	9.05
°C	50	60	70	80	90	
pH	8.98	8.93	8.90	8.88	8.84	