
Operating Instructions

Manual d'instruction

Model 8012
Dissolved Oxygen Probe

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The Model 8012 Dissolved Oxygen probe is designed to give a direct readout of dissolved oxygen concentrations in parts per million (ppm) in fresh water, when connected to a conventional pH/mV meter. Using a simple setting-up procedure the probe can be calibrated to measure concentrations in the range 0 to 14 ppm directly from the pH meter scale. The probe can be used for laboratory dissolved oxygen measurements, and is particularly suited for the B.O.D. (Biochemical Oxygen Demand) test, which estimates the amount of oxygen consumed by microscopic organisms in a water sample over a prescribed period.

The Model 8012 comprises an oxygen sensor capsule (a silver/lead galvanic cell) connected to a miniaturised circuit within the probe body. This circuit converts the current output of the capsule, produced by the electrochemical reduction of oxygen diffusing into it from the sample, to a millivolt signal compatible with the input requirements of a pH/mV meter. The circuit is powered from four tiny disposable cells contained in the probe body. Automatic temperature compensation is carried out using a thermistor mounted on the probe stem close to the sensor capsule. An adjustable control at the top of the probe body is used to calibrate the probe and meter combination.

The Model 8012 is supplied as part of a kit which includes two clip-on stirrer attachments, twelve cells, a spare sensor capsule, and a displacement funnel (with adaptor) for the B.O.D. bottle.

2.1 Requirements



In addition to the Model 8012 and the contents of the kit, the following will be required:

- A pH/mV with either a millivolt measuring range of at least 0 to 1400mV (i.e. so that 100mV corresponds to 1 ppm O₂) or a 0 to 14pH range (so that 1pH unit corresponds to 1 ppm O₂) with a back-off zero control capable of making the zero of the probe correspond to 0pH on the meter scale.
- A magnetic stirrer.
- A 5% solution of sodium sulphite (zero-concentration dissolved oxygen).
- A means of producing fully-aerated water at a controlled temperature.
- A mercury-in-glass thermometer, scaled 0 to 100°C or 0 to 50°C.

2.2 Contents of Kit

The complete contents of the Model 8012 kit (Fig 2.1) are:

| Qty | Description |
|-----|--|
| 1 | Model 8012 Probe body. |
| 2 | Stirrer bar attachments, (one with long spindle, one with short spindle). |
| 2 | Sensor capsules, in protective containers. |
| 1 | Displacement funnel (with adaptor) for the B.O.D. test - the funnel (and cover) may be used for B.O.D. bottles with narrow necks, whilst adding the knurled adaptor enables the funnel to be used on B.O.D. bottles with wide-tapered necks. |
| 12 | Disposable cells (in two sealed packs of six cells). |
| 1 | Instruction manual. |

2.3 Probe Assembly Instructions

Assembly of the Model 8012 is a simple operation comprising the two sequences described below, which should be followed carefully to give the best results.

2.3.1 Fitting the battery cells

(Fig. 2.3)

- a) Open the battery carrier compartment of the probe by drawing the transparent blue plastic cover down the stem of the probe, using a slight twisting action.
- b) Starting at the end of the carrier nearest the capsule, fit four cells into the cell carrier so that each cell has its positive pole (marked + on the cell) towards the top (cable) end of the probe. Each cell should snap neatly into position in the plastic carrier, and it should not be necessary to force any of the cells into place.

 **Warning.** Do not use any sharp pointed tools, e.g. tweezers, to remove the cells from the compartment: the cell casings may be punctured if handled carelessly.

- c) Slide the transparent blue plastic cover back over the compartment, making sure that the O-rings make a good seal.

 **Note.** Current drain from the cells is continuous as long as the cells are in the probe body, whether the probe is connected to a meter or not.


Fig. 2.3 Fitting battery cells in the probe

2.3.2 Fitting the Sensor Capsule

Caution. Clean and dry the area around the sensor capsule before the following instructions.


- a) Unscrew the old capsule from the sensor body and discard the capsule.
- b) Dry the sensor body with a paper tissue; ensure that the gold-coloured electrical contacts, and the recess into which the capsule screws, are clean and completely dry.


Caution.

- Ensure all the following steps are done carefully to avoid damaging the membrane covering the silver cathode.
- Do not leave the new capsule exposed to air for more than 30 minutes as the membrane will dry out.


- c) Access the new capsule by unscrewing the capsule housing


- d) Unscrew the capsule from the shorting cap.


- e) Remove and discard the rubber sealing ring from the new capsule.

- f) Dry the capsule with a paper tissue taking care not to damage the delicate, transparent membrane covering the silver cathode. Ensure that the gold-coloured contacts, and the threaded portion of the capsule are clean and completely dry.
- g) Fit the new 'O' ring supplied - see item e).
- h) Screw the capsule into the bottom end of the probe stem, making sure that the black rubber seal is compressed.
- i) Reverse the capsule housing and locate the lugs into the keyways on the capsule.
- j) Using the capsule housing tighten the the sensor as much as possible using finger pressure only to ensure a good seal.

- k) Clip onto the capsule, either the long or short spindle stirrer bar attachment. The choice of attachment depends upon the overall height of the sample bottle/ displacement funnel combination: for samples in a beaker, the short form is recommended.

3 MEASUREMENTS

3.1 Calibration of Probe

Before the probe can be used on samples to make measurements, it must first be calibrated by carrying out the following steps:

- a) Prepare in advance a solution of 5% sodium sulphite (5g/100ml of analytical grade sodium sulphite (Na_2SO_3) in demineralised water). Pour sufficient solution into a beaker to immerse the membrane end of the capsule fully when the stirrer attachment is fitted.
- b) Dip the probe into the solution, holding in position with a support stand. If this is the first zero check made on the probe after assembly, allow about 15 minutes for stabilisation to the zero oxygen level of the sulphite solution, otherwise a much shorter period should be sufficient, e.g. 4 to 5 minutes.
- c) Set the pH/mV meter to its 0-14pH scale and adjust the buffer/zero control so that the meter indicates exactly zero scale. This fixes the zero end of the meter scale to correspond to zero concentration oxygen. If there is insufficient adjustment on the control to permit this, as will be the case with most pH/mV meters (other than Models 7020, 7045 and 7046), the meter should be switched to the mV scale. It may be that on this scale adjusting the buffer/zero control will have no effect on the meter reading, but the small residual potential (a few mV) observed should have only a minimal effect on the overall accuracy of results obtained using the probe. This is explained in greater detail in the Appendix.
- d) Prepare a sample of fully-aerated water as follows:

Fill a beaker with demineralised water and aerate using either a small pump and aerator stone, or use an insufflator.

Whilst this method is sufficiently accurate for most applications, to obtain a precise 100% saturation solution, stir the beaker slowly, using a magnetic stirrer, so that the water is continually agitated but without breaking the liquid surface. Continue this stirring for approximately two hours to attain complete equilibrium. The water should be maintained at a temperature as close as possible to that of the sample (use the mercury thermometer to measure this). Note this temperature.
- e) Remove the probe from the sulphite solution and rinse thoroughly with demineralised water.
- f) Stand the beaker of aerated water on the stirrer and dip the probe into the water so that the end of the probe is immersed to a point above the thermistor housing on the side of the stem (this is to ensure that the thermistor can sense the water temperature accurately). Start the stirrer, and check that the stirrer bar is rotating at a rate which is as fast as possible whilst maintaining a smooth and steady rotation of the stirrer bar.
- g) Turn the scale length control knob on the top of the probe fully anti-clockwise and allow the pH/mV meter reading to stabilise for about five minutes.

* **Note.** It will be noted that if the probe is transferred to a solution containing oxygen after being immersed for some minutes in the sulphite solution (e.g. as occurs during the calibration sequence), then the response time of the probe will be somewhat longer than expected. This is quite normal, and for subsequent readings the probe will return to its move typical speed response of approximately 20s for 90% response to a step change in oxygen concentration. Where such a step-change is accompanied by a change in sample temperature slightly longer will be required for full stabilisation.

- h) When a steady reading is obtained, turn the scale length control on the probe slowly clockwise until the meter reading corresponds to the solubility value for the temperature noted in step d) above. This value is obtained from Table 1 in the Appendix at the end of this instruction manual. For the most accurate work a correction factor may be required to take into account variations in barometric pressure (e.g. an altitude correction). Table II in the Appendix lists correction factors against equivalent altitudes. As an example of how this may be calculated in practice, suppose the temperature of the aerated water is 10°C, and the barometric pressure is 750mm Hg. Then the control on the probe should be adjusted until the pH/mV meter indicates:

$$11.27 \times \frac{750}{760} = 11.12(\text{pH})$$

(or 1112mV on a
0-1400mV scale)

- i) The meter and probe combination is now calibrated to measure dissolved oxygen in the range 0 to 14ppm and is ready for measurement on samples.

* **Note.** When not in use the probe should always be stored so that the sensor capsule is in water. The capsule cannot be reactivated once it has been allowed to dry out.

3.2 Measurements on Samples

To carry out measurements on samples perform the following steps:

- Stand the beaker of sample on the magnetic stirrer and immerse the probe in the sample, ensuring that the thermistor housing is completely covered.
- Switch on the stirrer using the same stirring rate as for the aerated solution, and check that the stirrer bar rotates smoothly. Ensure that no air bubbles are trapped around the membrane of the sensor capsule.
- Allow sufficient time for the probe to stabilise to the sample temperature, and to provide a steady reading on the pH/mV meter.
- Read off the dissolved oxygen concentration directly from the meter scale.

* **Note.** Since the output from the probe depends upon oxygen diffusing across the gas-permeable membrane, it is essential that the aqueous layers on the outside of the membrane do not become depleted of oxygen. It is thus clear that efficient stirring of the liquid in the vicinity of the membrane is necessary, and this requirement is met by the stirrer bar attachment of the probe. Two attachments are provided, with either short or long spindles, depending on the particular type of flask or beaker used to hold the sample. Since there is some dependence on output on the stirring rate, it is recommended that the rate used to calibrate the probe is also used for all measurements.

- e) After completing a measurement or series of measurements, store the probe so that the sensor capsule is in water, as recommended in section 4.

3.3 B.O.D. Measurement (Fig. 3.1)

To facilitate measurements in the B.O.D. (Biochemical Oxygen Demand) test, the Model 8012 is supplied with a specially-designed displacement funnel, adaptor and cover for use with most sizes and neck tapers of B.O.D. bottle. According to the type of bottle used, either the displacement funnel and cover alone are used to seal the bottle (narrow taper), or the funnel, cover and knurled adaptor (wide taper).

In both cases the B.O.D. bottle is filled to the brim before inserting the funnel. The probe is then inserted through the cover, which holds the probe securely in position whilst measurement is being made. The displacement funnel retains the water displaced as the probe is immersed in the sample. The displaced water flows back into the bottle when the probe is removed. The bottle thus remains completely full, and since the probe consumes a negligible amount of oxygen, initial and final measurements for the B.O.D. test can be made on the same bottle of sample. For example, at 20°C the probe consumes approximately 0.5% of the available oxygen in a closed 250ml sample in one hour. A measurement, however, takes only a few minutes, so that the consumption is negligible.

It should be noted that, in comparison with the Winkler titration method, which ensures sterilisation of the bottles automatically by the presence of iodine, the probe measurement technique will require sterilisation of the B.O.D. bottles after each test, before they can be reused. A procedure for sterilising B.O.D. bottles is as follows:

Prepare 1 litre of 1% m/v sulphuric acid and add 2.5g iodine and 12.5g potassium iodide. Shake to dissolve. (This solution should be discarded when the brown colour fades.) The bottles should be rinsed with 5 to 10ml of this solution, shaking well to coat the bottle walls. Stand for 15 minutes, pour away the wash solution and rinse thoroughly with water and finally with distilled or de-ionised water.

3.4 Storage

When the probe is not being used for measurements, it should be stored so that the membrane is in water. It should never be allowed to dry out – the sensor capsule cannot be reactivated if this happens.

Remember that as long as the cells are inserted in the cell carrier compartment, there will be a current drain from them, whether the probe is connected to a pH/mV meter or not.

If the probe is not required for an extended period, i.e. longer than one month, it is recommended that the cells and sensor capsule be removed and stored separately. To remove the cells from the probe, first slide down the transparent plastic cover of the carrier compartment and press the flat bar marked 'EJECTOR' on the back of the compartment.

Warning. Do not use any sharp pointed tools, e.g. tweezers, to remove the cells from the compartment: the cell casings may be punctured if handled carelessly.

Always take care on replacing cells in the probe to ensure correct polarity - the positive pole (marked + on the cell) should be towards the top (cable) end of the probe (see Fig. 2.3).

After unscrewing the sensor capsule from the probe, the stirrer attachment is removed by inserting a screwdriver blade into the clip-on fitting and gently prising it away from the sensor capsule end.

Replace the shorting cap provided on the capsule, retaining the sealing washer, and return the capsule to its protective vial. Ensure that the sponge rubber pad in the vial is kept moistened with water.


Fig. 3.1 Using the Model 8012 with B.O.D. Bottles

4 ROUTINE MAINTENANCE & STORAGE

4.1 Replacement of Sensor Capsule and Cells

If after a period of extended operation the probe becomes very sluggish and slow to respond to changes in oxygen concentration, or cannot be calibrated, then the capsule should be replaced. Difficulties experienced in calibration could also indicate that the cells voltage is low: eject the cells and replace all of them with a new set. The cells supplied with a probe are a Ray-O-Vac silver oxide type RS 13-G or RW 48-X (Part No. 0231 138), and these are the recommended cells to be used as spares. They are available in packs of six cells from ABB Kent-Taylor, or suitable equivalent mercury type or silver-oxide types may be obtained from stockists of hearing aid/digital watch cells. The following is a list of some possible equivalent cells:

| Manufacturer | Mercury type | Silver oxide type |
|--------------------|--------------|--------------------|
| Ray-O-Vac | R13 | RS 13-G or RW 48-X |
| Mallory | RM 13H | D393 |
| Berec (Ever-Ready) | RM 13H | B/SR 48/H |
| Varta | V13 HM | V13 HS |

The expected life of these cells in continuous operation (i.e. inserted in the probe) is about three months: the mercury types have slightly longer operational life, whilst the silver oxide types have a longer shelf life.

Warning. NONE of these cells are rechargeable types, and recharging of the cells must **NEVER** be attempted under any circumstances.

5 SPARES & ACCESSORIES

| Title | Part No. |
|--|----------|
| Silver oxide cell Ray-O-Vac type RS 13-G or RW 48-X (pack of six) | 0231 138 |
| Sensor capsule, complete with shorting cap, sealing washer, in protective vial | 8012 170 |
| Shorting Cap and Tightening Tool | 8012 650 |
| Sealing washer for sensor capsule | 0211 025 |
| Stirrer bar attachment, long spindle | 8012 130 |
| Stirrer bar attachment, short spindle | 8012 120 |
| Displacement funnel | 8012 530 |
| Displacement funnel cover | 8012 540 |
| Funnel adaptor (U.K. version) | 8012 550 |
| Funnel adaptor (U.S. version) | 8012 560 |


6.1 Sensor capsule

The Model 8012 Dissolved Oxygen Probe depends for its operation on a silver/lead galvanic cell (the sensor capsule), in which an electrolyte is retained by a gas-permeable membrane. Oxygen which diffuses through this membrane is reduced at the silver cathode of the capsule to hydroxide ions, the necessary electrons being supplied from a silver surface: at the same time hydroxide ions in the neighbourhood of the lead anode react with it, giving up their charge, and thus a current is caused to flow in a circuit connected to the probe.

At the silver cathode, the reaction can be represented by the equation:

and at the lead anode, by the equation:

The magnitude of the current is directly proportional to the diffusion rate of oxygen across the membrane, and this is in turn proportional to the partial pressure of oxygen in the sample and the permeability of the membrane. Unfortunately, the permeability is not constant, but varies exponentially with temperature; thus it is necessary to correct for this effect. This can be conveniently achieved using an appropriate thermistor. In addition, the thermistor is used to compensate for the solubility variation of oxygen with temperature, and this enables concentration measurement to be directly read off the pH meter scale.

Since the lead anode is consumed by the reaction, the sensor capsule has a finite life, depending upon the amount of lead initially available, and the capsule has been designed to both minimise oxygen consumption (particularly important in the B.O.D. test, for example) and to provide as much lead as possible, to obtain the longest possible operating life for each capsule.

6.2 Electronic circuit

The electronic circuit, powered by the four disposable cells, incorporates a micropower operational amplifier, and serves to convert the current output of the sensor capsule to a millivolt signal compatible with pH meter requirements. A miniature encapsulated bead thermistor is fitted to the probe near the sensor capsule, and provides temperature compensation by sensing sample temperature variations and altering the sensitivity of the amplifier accordingly.

The amplifier components are mounted on a narrow miniature circuit board located in the centre of the probe stem, except for a 'scale length' control (connected by flexible circuitry), positioned at the top of the probe. To save space, miniaturised components are used mounted on the board, with a filter network to separate it from the amplifier stage.

The disposable cells provide a positive and negative supply to the amplifier, and this occurs continuously as long as the cells are inserted in the probe.

Measurement Range:

0 to 14ppm (mg/l-1) dissolved oxygen.

Temperature Range:

0 to 40°C.

Accuracy:

± 0.2ppm at ± 5°C of initial calibration temperature.

Repeatability:

± 0.1ppm.

Response Time:

Typically 20 seconds for 90% of a step-change of oxygen concentration at 20°C, slightly longer when a sample temperature also varies.

Temperature Compensation

Automatic correction by means of integral thermistor in probe.

Cells:

Four cells, (nominal cell voltage 1.2V).

Dimensions:

Probe – 184mm (7.25in.) long approximately excluding stirrer bar attachment, 12.7mm (0.5in.) diameter.

Kit – 250 x 150 x 65mm approximately overall (9.88 x 6 x 2.56 in).

Weight of Kit:

0.8kg (1.25lb).

A.1 Solubility of Oxygen in Pure Water

Table I below gives the solubility of oxygen in ppm (mg/l-1) in pure water at equilibrium with water vapour saturated air at the standard atmosphere pressure of 760mm mercury.

Note. Table A.1 is abstracted

from Table IVb of 'International Oceanographic Tables' Vol. 2, National Institute of Oceanography of Great Britain and UNESCO 1973 (0-34°C) and R.F. Weiss, Deep-Sea Res., 1970, 17, 721, (36-40°C).

For different barometric pressures (e.g. at different altitudes) a correction should be applied by multiplying the solubility figure given above for the given sample temperature by the appropriate factor. See the Table A.2.

Assuming an atmospheric pressure at sea level of 760mm Hg, the correction factors for altitude variations shown in Table A.2 should be applied to the solubility figures given in Table A.1.

For the most accurate results it is recommended to measure the local barometric pressure and calculate the correction factor from the fraction $P/760$, where P is the barometric pressure in mm Hg. e.g. for a barometric pressure of 740mm Hg, the solubility figures in Table A.1 should be multiplied by $740/760 = 0.974$.

A.2 Effects of Single Point Calibration Method on Accuracy of Measurement

As previously stated in Section 3.1 (Step c) it is possible when using some types of pH/mV meter that the instrument when used on the mV scale cannot be precisely set to zero when the probe is immersed in sodium sulphite solution, because the buffer/zero control on the meter has no effect on the meter reading. In this event zero concentration oxygen in the sodium sulphite solution will produce a small reading on the meter of a few millivolts instead of zero scale (0mV).

Fig. A.1 illustrates this possibility for a sample where the calibration of the meter is made at 20°C and 760mm Hg. From Table A.1 the solubility is read off as 9.07 ppm O₂, and a probe/meter combination in which calibration is possible with full

back-off adjustment of the buffer/zero control is represented by Line A on the graph. Line B shows a calibration where the residual voltage in a zero ppm O₂ solution was not backed-off by the buffer/zero control of the meter. It can be seen that the two lines diverge (and thus the error increases) away from the calibration point.


However, this is only apparent at measurement significantly removed from the calibration point i.e. at very low oxygen concentrations, and at no point does it exceed the overall accuracy of the probe, as stated in the Section 7. For most requirements, the small residual potential (in this example 20 mV approximately) can be ignored without detriment to the overall accuracy of measurement.

| Temp °C | ppm O ₂ | Temp °C | ppm O ₂ |
|---------|--------------------|---------|--------------------|
| 0 | 14.59 | 21 | 8.90 |
| 1 | 14.19 | 22 | 8.73 |
| 2 | 13.81 | 23 | 8.55 |
| 3 | 13.44 | 24 | 8.40 |
| 4 | 13.08 | 25 | 8.24 |
| 5 | 12.75 | 26 | 8.08 |
| 6 | 12.42 | 27 | 7.94 |
| 7 | 12.12 | 28 | 7.80 |
| 8 | 11.82 | 29 | 7.66 |
| 9 | 11.54 | 30 | 7.54 |
| 10 | 11.27 | 31 | 7.41 |
| 11 | 11.01 | 32 | 7.28 |
| 12 | 10.75 | 33 | 7.15 |
| 13 | 10.52 | 34 | 7.04 |
| 14 | 10.28 | 35 | 6.93 |
| 15 | 10.07 | 36 | 6.82 |
| 16 | 9.85 | 37 | 6.71 |
| 17 | 9.64 | 38 | 6.61 |
| 18 | 9.44 | 39 | 6.51 |
| 19 | 9.25 | 40 | 6.41 |
| 20 | 9.07 | | |

Table A.1 Solubility of Oxygen in Pure Water

| Approximate Altitude | | Correction Factor |
|----------------------|-----------|-------------------|
| Metres | Feet | |
| | Sea-level | 1.000 |
| | 500 | 0.982 |
| | 1000 | 0.965 |
| | 1500 | 0.948 |
| 500 | | 0.942 |
| | 2000 | 0.931 |
| | 2500 | 0.915 |
| | 3000 | 0.899 |
| 1000 | | 0.887 |
| | 3500 | 0.883 |
| | 4000 | 0.867 |
| | 4500 | 0.851 |
| | 5000 | 0.836 |
| 1500 | | 0.835 |
| | 5500 | 0.821 |
| | 6000 | 0.806 |
| | 6500 | 0.791 |
| 2000 | | 0.784 |
| | 7000 | 0.777 |
| | 7500 | 0.762 |
| | 8000 | 0.748 |
| 2500 | | 0.737 |
| | 8500 | 0.735 |
| | 9000 | 0.721 |
| | 9500 | 0.708 |
| | 10000 | 0.694 |
| 3000 | | 0.692 |

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
REMARQUES

REMARQUES

ABB has Sales & Customer Support expertise
in over 100 countries worldwide

www.abb.com

The Company's policy is one of continuous product
improvement and the right is reserved to modify the
information contained herein without notice.

© ABB 2002

Printed in the UK (08.02)

The ABB logo consists of the letters 'A', 'B', and 'B' in a bold, black, sans-serif font. The 'A' is slightly larger than the 'B's, and they are all closely spaced together.

ABB Limited

Howard Road, St. Neots
Cambridgeshire, PE19 8EU
UK

Tel: +44 (0)1480 475 321

Fax: +44 (0)1480 217 948

ABB Inc

2175 Lockheed Way
Carson City, NV 89706
USA

Tel: +1 775 883 4366

Fax: +1 775 883 4373

ABB Automation

100 Rue de Paris
F-91342 Massy Cedex
France

Tel: +33 1 64 47 20 00

Fax: +33 1 64 47 20 16