Operating Instructions

Manual d'instruction

Model 8001
Fluoride Electrode

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The Company

We are an established world force in the design and manufacture of instrumentation for industrial process control, flow measurement, gas and liquid analysis and environmental applications.

As a part of ABB, a world leader in process automation technology, we offer customers application expertise, service and support worldwide.

We are committed to teamwork, high quality manufacturing, advanced technology and unrivalled service and support.

The quality, accuracy and performance of the Company's products result from over 100 years experience, combined with a continuous program of innovative design and development to incorporate the latest technology.

The NAMAS Calibration Laboratory (No. 0255) is just one of ten flow calibration plants operated by the Company, and is indicative of our dedication to quality and accuracy.

EN ISO 9001:2000

Cert. No. Q5907

EN 29001 (ISO 9001)

Lenno, Italy – Cert. No. 9/90A

Stonehouse, U.K.

Health and Safety

To ensure that our products are safe and without risk to health, the following points must be noted:

1. The relevant sections of these instructions must be read carefully before proceeding.
2. Warning labels on containers and packages must be observed.
3. Installation, operation, maintenance and servicing must only be carried out by suitably trained personnel and in accordance with the information given.
4. Normal safety precautions must be taken to avoid the possibility of an accident occurring when operating in conditions of high pressure and/or temperature.
5. Chemicals must be stored away from heat, protected from temperature extremes and powders kept dry. Normal safe handling procedures must be used.
6. When disposing of chemicals ensure that no two chemicals are mixed.

Safety advice concerning the use of the equipment described in this manual or any relevant hazard data sheets (where applicable) may be obtained from the Company address on the back cover, together with servicing and spares information.

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1 INTRODUCTION

The Model 8001 fluoride electrode senses activity and hence concentration of fluoride ion in aqueous solutions from saturated fluoride solutions to 10^{-6}M (approx. $0.02\text{mg l}^{-1} \text{ F}^{-}$). The electrode should be used with a calomel reference electrode (Part No. 1431-210 or 1370-210) and an expanded scale pH meter such as Model 7045/6. When used with a Model 7046 with the function switch set to monovalent anion, the concentration of fluoride ions in the sample may be read directly from the meter scale, thus eliminating the need for a calibration graph.

The fluoride electrode has as a sensing tip a single crystal of lanthanum fluoride which is mounted in a highly chemically resistant polyolefin body. Electrical contact is made with an internal electrolyte which is permanently sealed inside the body.


Fig. 1.1 Model 8001 Fluoride Electrodes

2 GENERAL INFORMATION

2.1 Requirements

- Model 8001 fluoride ion-selective electrode.
- An expanded scale pH meter (such as Model 7046, which has direct reading concentration scales).
- A calomel reference electrode (Part No. 1431-210 or 1370-210).
- Standard fluoride solutions (see Section 3).
- Ionic strength and pH adjustment buffer (see Section 3).
- Magnetic stirrer and stirrer bars.
- Plastic labware.

2.2 Accessory Equipment (Electrode Holder, Part No. 2610 300)

This item is not supplied with the kit. It may be used on a laboratory stand (Part No. 2610 130) to support both the fluoride and reference electrodes.

2.3 Storage

Model 8001 fluoride electrode can be immersed in a dilute fluoride solution when not in use. For long periods the electrode should be stored dry with the protective end cap replaced to prevent the crystal being scratched or damaged. Before use the electrode tip should be rinsed in distilled or de-ionised water.

3 MEASUREMENT TECHNIQUES

3.1 Buffer Solution

The pH of the sample and standard solutions should be above pH 5 to prevent complexation of fluoride with hydrogen ions and below pH 7 to avoid hydroxide interference (see Section 4.3).

The majority of fluoride measurement applications (e.g. potable water) require use of a Total Ionic Strength Adjustment Buffer (TISAB).

Addition of TISAB ensures constant ionic strength between samples and standards and adjusts the pH to 5.0 - 5.5.

TISAB will also complex iron and aluminium preferentially, breaking up fluoride complexes with these metal ions.

In a 1mg l^{-1} fluoride sample, TISAB complexes about 5mg l^{-1} aluminium or iron.

To prepare 1 litre of TISAB:

- a) Place 500ml of distilled or de-ionised water in a 1 litre beaker. Add $57(\pm 1)\text{ml}$ of glacial acetic acid, $58(\pm 0.5)\text{g}$ sodium chloride and $4(\pm 0.5)\text{g}$ CDTA (trans-1, 2-diaminocyclohexane-N,N',N'-tetra-acetic acid). Stir to dissolve. Place the beaker in a water bath and add carefully sufficient 5M NaOH (200g l^{-1}) to adjust the pH to between 5.0 and 5.5, using a pH electrode pair to monitor the pH of the mixture. Transfer to a one litre volumetric flask and dilute to the mark with distilled or de-ionised water.

TISAB is effective with a wide variety of different samples of different ionic strengths when used at a sample to buffer ratio of 1:1. However, in many instances (e.g. most potable waters in the U.K.), it will be found effective when used at a sample to buffer ratio of 10:1.

If the concentration of aluminium or iron in the sample is appreciable, i.e. greater than about 5mg l^{-1} , then an alternative buffer solution should be used. To prepare 1 litre proceed as follows:

- b) Place 500ml of distilled or de-ionised water in a 1 litre beaker. Add carefully $84 (\pm 1)\text{ml}$ of concentrated hydrochloric acid (1.18s.g.) and mix well. Dissolve in this solution $242(\pm 1)\text{g}$ tris (hydroxymethyl) methylamine and $230 (\pm 1)\text{g}$ sodium tartrate ($[\text{CHOH.COO Na}]_2 \cdot 2\text{H}_2\text{O}$). Ensure that the solution is at room temperature and then transfer to a one litre volumetric flask. Dilute to the mark with distilled or de-ionised water.

Warning. Appropriate safety precautions must be taken when handling concentrated hydrochloric acid.

When used at a sample to buffer ratio of 1:1 this buffer solution complexes up to 100mg l^{-1} aluminium and iron at a fluoride concentration of 1mg l^{-1} . For higher levels of aluminium and iron the more complicated buffering procedure in reference 1 may be useful.

3.2 Standard Solutions

Warning. Sodium and potassium fluoride solutions are poisonous, and must always be handled with great care.

The standard solutions for calibration of the electrode can conveniently be prepared by serial dilution of fresh standard fluoride solutions, e.g.

0.1M sodium fluoride:
 $4.199(\pm 0.001)\text{g}^{-1}$
analytical reagent grade NaF

or

0.1M potassium:
 $5.810 (\pm 0.001)\text{g}^{-1}$
fluoride analytical reagent grade KF

or
1000 mg l^{-1} fluoride ions:
 $2.210 (\pm 0.001)\text{g}^{-1}$
analytical reagent grade NaF
or
1000 mg l^{-1} fluoride ions
 $3.058 (\pm 0.001)\text{g l}^{-1}$
analytical reagent grade KF

These standard solutions should always be stored in plastics labware, as fluoride reacts with glass. Before measurement dilute standards should be treated with buffer as previously described.

Low concentration fluoride solutions, e.g. less than 10^{-4}M (1.9mg l^{-1}) should be prepared fresh daily.

3.3 General Procedure

- 1) Rinse the reference and fluoride electrodes thoroughly in de-ionised water. Wipe dry with a tissue.
2. Insert reference and fluoride electrode connectors into the appropriate sockets on the meter.
- 3) Place both electrodes in appropriate standardising solutions and read the potential obtained after the reading has stabilised. Rinse and carefully wipe the electrodes with a tissue between solutions to prevent carryover. Stirring of both standards and samples is recommended, preferably using a magnetic stirrer and stirrer bar.
- 4) Prepare a calibration graph by plotting potentials on the linear axis or linear/logarithmic graph paper (see Fig. 3.1).
- 5) Place both electrodes in the unknown solution and read the potential. Sample and standards should be at the same temperature. Interpolate sample activity or concentration from the calibration graph.

3.4 Notes and Recommendations

- 1) It is strongly recommended that all measurements be carried out in plastics labware. Use of glass may result in erratic or slow responses.
- 2) If using Model 7046 pH meter and measurements fall within the Nernstian range of the electrode, calibrate the concentration scale directly and dispense with the calibration curve.
- 3) For accurate results the composition of the standards should be as similar as possible to that of the samples. Thus when the composition of the samples is well defined, standards should be prepared from model samples with known amounts of fluoride ions.

- 4) In common with many other ion-selective electrodes the fluoride electrode reaches equilibrium more quickly after an increase in determinand concentration than after a decrease.

Accordingly, it is recommended that the electrode is calibrated in the standard solutions in ascending order of concentration.

3.5 The Known Addition Method

There are some advantages of using the Known Addition Method of analysis for occasional samples with fluoride concentrations greater than the lower Nernstian limit in the experimental conditions used.


Fig. 3.1 Model 8001 Fluoride Electrode Calibration Graph

This reduces the number of standard solutions required to one. The concentration of this standard solution should be between ten and a hundred times that of the sample. The experimental slope factor (S_{mV}) of the electrode should be known.

To analyse a sample by this method, pipette an aliquot, V_s ml. of sample into a beaker and measure the potential in this solution E_1 mV. Add a known volume, V_a ml, of the standard (where $10V_a \approx V_s$) of concentration C_a mol l^{-1} and after stirring measure the new potential of the electrodes, E_2 mV. It should be arranged that the potential change is not less than 5mV.

The concentration of fluoride in the original sample, C_s mol l^{-1} , may then be calculated from the formula:

$$C_s = \frac{C_a V_a (V_a + V_s)^{-1}}{\text{antilog}_{10}[(E_2 - E_1)S^{-1}] - V_s (V_s + V_a)^{-1}}$$

3.6 Electrode Malfunction

Caution. Do not disassemble a suspected faulty electrode – send back to for examination.

- 1) Noisy or unstable response:
 - a) Check connections.
 - b) Check for air bubbles trapped on crystal surface.
 - c) Check for surface deposits on crystal.
 - d) Check for reference electrode liquid junction partially or completely blocked.
 - e) Check whether samples and standards are buffered with TISAB.
- 2) Electrode acts as open circuit (readout on instrument stays off-scale or swings erratically about):
 - a) As 1a), b), c).
- 3) Slow or no response to changes in activity/concentration:
 - a) As 1c).
 - b) Samples not adequately buffered.
 - c) Solutions contain high concentrations of carboxylate anions (see Section 4.1).
- 4) Readings drifting:
 - a) Membrane dirty – wipe with methanol or lightly polish with fine (e.g. 0.3μ particle size) alumina powder.
 - b) Changes in sample temperature – ensure stable temperature.
 - c) Faulty reference electrode. Change to a new one.
 - d) Glassware used – use plastics labware.
 - e) As 3 c).
- 5) Low slope on calibration curve:
 - a) Check that standard solutions are fresh and have not been contaminated.
 - b) As 1 c), 4 a).
 - (c) As 4 b).

4 ELECTRODE CHARACTERISTICS

4.1 Response range

Fluoride concentrations down to 10^{-6}M F^- (0.019mg l^{-1}) can be measured in neutral solutions and in the absence of complexing agents provided adequate precautions are taken to prevent contamination of dilute solutions (see Section 3.2). The upper limit of detection is a saturated fluoride solution.

4.2 Temperature Effects

Since the standard potentials of both fluoride and reference electrodes are affected by temperature, standards and samples should be at the same temperature to within $\pm 1^\circ\text{C}$. The slope of the fluoride electrode also varies with temperature according to the factor $2.3RT/F$ in the Nernst equation.

The fluoride electrode can be used at temperatures from 0 to 60°C provided that thermal equilibrium has occurred. Where substantial temperature change takes place, up to one hour should be allowed to re-establish equilibrium.

4.3 pH effects

Below pH5 a portion of the fluoride is complexed with hydrogen ions to form HF and HF_2^- , neither of which can be detected by the fluoride electrode. Above pH7 hydroxide ions interfere. TISAB buffers the solution to between 5.0 and 5.5 pH to avoid complexation of fluoride and hydroxide interference (for preparation of TISAB, see Section 3).


Fig. 4.1 Model 8001 Fluoride Electrode pH Response in Standard Solutions

4.4 Complexation

Aluminium, iron [Fe(III)], calcium, magnesium and other polyvalent cations, as well as hydrogen, complex with fluoride. The degree of complexation depends on various factors such as the pH and total ionic strength of solution and the concentration of complexant and fluoride.

CDTA which is present in TISAB reagent preferentially complexes these polyvalent cations in the sample – see Section 3.

4.5 Interferents

The major interferent in the determination of fluoride is hydroxide ion. In this case the electrode responds to hydroxide ion in solution although the response is not Nernstian. This effect can be eliminated by buffering below pH7.

4.6 Response time

Stirring of the sample is strongly recommended to accelerate response time. Typical responses vary, from several seconds in 10^{-1}M to 10^{-4}M fluoride, to about one to two minutes in 10^{-6}M solution where the electrode is near its limit of detection. Response times tend to be slightly quicker for increases in concentration than for decreases.

At low fluoride concentrations the presence of large quantities of carboxylate anions result in sluggish response times (see reference 2).

The fluoride electrode consists of a single crystal of lanthanum fluoride which is bonded in an inert polyolefin body. The crystal is an ionic conductor in which only fluoride ions are mobile. Inside the electrode an aqueous filling solution, containing fluoride and chloride ions, ensures a constant potential between the membrane and the internal Ag/AgCl reference element.

When the electrode is in contact with a fluoride solution, a potential, which is measured against an external constant reference potential, develops across the membrane. Changes in this potential are caused by changes in fluoride activity in the sample.

The Nernst equation describes the potential behaviour of the electrode in pure fluoride solutions:

$$E = E_o - \frac{2.3RT}{F} \log_{10} a_{F^-}$$

where

E = measured potential of system

E_o = conditional standard potential, which varies according to the choice of reference electrode

$2.3RT/F$ = Nernst factor (59.2mV at 25°C) R and F are constants and T is the temperature in Kelvin

a_{F^-} = fluoride activity in the sample solution

The activity of fluoride ions, a_{F^-} , is related to their concentration $[F^-]$, by the fluoride ion activity coefficient, γ_{F^-} as follows:

$$a_{F^-} = \gamma_{F^-} [F^-]$$

This activity coefficient will remain virtually constant if the ionic strength is held constant. In practice this is usually achieved by adding a constant level of background electrolyte to all standards and samples.

Thus if γ_{F^-} is a constant

$$E = E_o - \frac{2.3RT}{F} \log_{10} [F^-]$$

Interconversion of concentration units

$$\begin{aligned} 1.00 \times 10^{-3} M F^- &= 1.00 \times 10^{-3} \text{mol l}^{-1} \\ &= 19.0 \text{mg l}^{-1} \\ 1.00 \text{mg l}^{-1} F^- &= 1.00 \text{ppm} \\ &= 5.26 \times 10^{-5} \text{mol l}^{-1} \end{aligned}$$

6 SPECIFICATION

Concentration range:

Saturated solutions to $10^{-6} M F^-$
(0.019mg l⁻¹)

Lower Nernstian Concentration Limit:

$5 \times 10^{-6} M F^-$

Temperature range:

0 to 60°C (constant use)

60 to 70°C (intermittent use)

Electrode Resistance – typically 100K:

pH range (see Fig. 4.1)

pH 5 to 7 at $10^{-6} M$ (0.019mg l⁻¹ F⁻)

pH 5 to 11 at $10^{-1} M$ (1900mg l⁻¹ F⁻)

Reproducibility:

Better than $\pm 2\%$ of concentration

Minimum Sample Size

5ml in 50ml beaker

Storage

Store dry or in dilute fluoride solution (see Section 3.2). For long-term storage replace protective cover.

Dimensions

	White Cap	Red Cap
Lead length:	1m	1m
Overall length:	138mm	165mm
Stem diameter:	12mm	12mm
Length under cap:	105mm	120mm

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Products

Automation Systems

- *for the following industries:*
 - Chemical & Pharmaceutical
 - Food & Beverage
 - Manufacturing
 - Metals and Minerals
 - Oil, Gas & Petrochemical
 - Pulp and Paper

Drives and Motors

- *AC and DC Drives, AC and DC Machines, AC motors to 1kV*
- *Drive systems*
- *Force Measurement*
- *Servo Drives*

Controllers & Recorders

- *Single and Multi-loop Controllers*
- *Circular Chart and Strip Chart Recorders*
- *Paperless Recorders*
- *Process Indicators*

Flexible Automation

- *Industrial Robots and Robot Systems*

Flow Measurement

- *Electromagnetic Flowmeters*
- *Mass Flow Meters*
- *Turbine Flowmeters*
- *Flow Elements*

Marine Systems & Turbochargers

- *Electrical Systems*
- *Marine Equipment*
- *Offshore Retrofit and Refurbishment*

Process Analytics

- *Process Gas Analysis*
- *Systems Integration*

Transmitters

- *Pressure*
- *Temperature*
- *Level*
- *Interface Modules*

Valves, Actuators and Positioners

- *Control Valves*
- *Actuators*
- *Positioners*

Water, Gas & Industrial Analytics Instrumentation

- *pH, conductivity, and dissolved oxygen transmitters and sensors*
- *ammonia, nitrate, phosphate, silica, sodium, chloride, fluoride, dissolved oxygen and hydrazine analyzers.*
- *Zirconia oxygen analyzers, katharometers, hydrogen purity and purge-gas monitors, thermal conductivity.*

Customer Support

We provide a comprehensive after sales service via our Worldwide Service Organization. Contact one of the following offices for details on your nearest Service and Repair Centre.

United Kingdom

ABB Limited
Tel: +44 (0)1453 826661
Fax: +44 (0)1453 829671

United States of America

ABB Inc.
Tel: +1 775 850 4800
Fax: +1 775 850 4808

Client Warranty

Prior to installation, the equipment referred to in this manual must be stored in a clean, dry environment, in accordance with the Company's published specification. Periodic checks must be made on the equipment's condition.

In the event of a failure under warranty, the following documentation must be provided as substantiation:

1. A listing evidencing process operation and alarm logs at time of failure.
2. Copies of all storage, installation, operating and maintenance records relating to the alleged faulty unit.

ABB has Sales & Customer Support
expertise in over 100 countries worldwide

www.abb.com

The Company's policy is one of continuous product
improvement and the right is reserved to modify the
information contained herein without notice.

Printed in UK (11.04)

© ABB 2004

ABB Limited

Oldends Lane Stonehouse
Gloucestershire
GL10 3TA
UK

Tel: +44 (0)1453 826661
Fax: +44 (0)1453 829671

ABB Inc.

Analytical Instruments
9716 S. Virginia St., Ste. E
Reno, Nevada 89521
USA

Tel: +1 775 850 4800
Fax: +1 775 850 4808